

**ViggoWAY™**  
**3 Way Stop Cock**

**MATERIALS USED:**  
 •Polycarbonate (PC), HDPE, PP

**COATING MATERIAL :**  
 •Silicone Dispersion

**INDICATIONS:**  
 •To deliver simultaneous infusion of two different IV fluids  
 •To administer bolus IV injections through existing access site for IV infusion  
 •To take blood samples through existing IV access site, without a second venipuncture  
 •For needle-free access in high-risk patients

**CONTRAINDICATIONS:**  
 •Product should not be used in patients with known hypersensitivity to any of the materials used including coating materials  
 •Product should not be used for Administration of highly viscous fluids  
 •Product should not be used for blood transfusion














**INSTRUCTIONS FOR USE:**  
 •The connecting port of this device is in compliance to EN 20594-1 & ISO 594-2  
 •Attachment of three-way stop cock to IV access  
 •Ensure that venous access (e.g. IV cannula) is patent by flushing with saline  
 •Inspect packing of 3-way Stop Cock visually to ascertain that package is intact  
 •Remove the Stop Cock from sterile packing taking care to maintain sterility  
 •Connect one or two IV Infusion lines to the ports as required  
 •Remove protector from male portion of stop cock and hold it upright  
 •Open clamps of IV infusion lines and displace air entirely by rotating the tap to open position (all 3 flow directions should be open)  
 •Rotate the tap to closed position (output from stop cock is closed)  
 •Connect it to venous access device (e.g. IV cannula) and secure the connection with rotating transparent lock  
 •Rotate the tap to open position and adjust the flow of infusion fluid

**REMOVAL OF 3-WAY STOP-COCK**  
 •Rotate the tap to closed position such that connection with IV access is terminated.  
 •Clamp input IV infusion lines  
 •Keep injection stopper/other device ready for closing IV access (e.g. IV cannula) after removal  
 •Unscrew the transparent rotator lock and remove 3-way stop cock; close IV access with injection stopper or other device

**TERM OF USE :**  
 •Recommended Maximum Duration of Use: Depends on type of medication and other factors; refer to your facility's guidelines or CDC Guidelines

**WARNINGS:**  
 •The use of this product is restricted to a qualified doctor or a paramedic  
 •Read instructions before use  
 •The product should be used according to the instructions for use  
 •If there is any change in expected performance of the device or in case of any malfunction the device should be immediately removed & sent back to supplier for analysis  
 •For known/reported adverse events associated with use of this device, refer to the Clinical Evaluation Report HH-QA-CER-3WSC  
**VIGGO MEDICAL DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE**  
 •Any device that is connected to this product must comply with EN 20594-1 & ISO 594-2 in order to achieve the intended performance of this product & to avoid leakage in this connection  
 •The product should not be reprocessed  
 •Visually inspect and carefully check the product and packaging before use. Improper transport and handling may cause structural and/or functional damage to device or packaging  
 •The product is guaranteed sterile & non-pyrogenic if the package has not been opened or damaged  
 •Do not clean or re-sterilise. FOR SINGLE USE ONLY. Discard after use  
 •Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or infections  
**Disposal/Discard:** Dispose off/Discard the used Device in accordance with your Country's Healthcare and Safety Regulations

**INSTRUCTIONS FOR USE**

	Product Ref. No.
	Batch Number
	Date Of Manufacturing
	Use By
	Do Not Reuse
	Sterilised by Ethylene Oxide Gas
	Caution, Consult Accompanying Documents
	Consult Instructions for Use
	Non-Pyrogenic
	Do Not Use if Package is Damaged
	Do not re-sterilize
	Keep away from sunlight
	Keep Dry

