

An Integrated PCR-Microarray Lab-on-Chip Assay for Detection and Identification of 4 Category A Biological Agents (Anthrax, Plague, Small Pox, Tularemia)

Kian-Leong Ong¹, Monica Capozzoli², Gregory Goh¹, Audrey Fong¹, Cheryl Wang¹, Rosemary Tan¹

Veredus Laboratories Pte Ltd, 83 Science Park Drive, #04-02, The Curie, Singapore Science Park, Singapore 118258¹; STMicroelectronics, Molecular Diagnostics R&D, Microfluidics Division, Computer Peripherals Group, Stradale Primosole 50, 95121 Catania, Italy²

Introduction

Current fears over the threat of bioterrorism attacks using Category A biological agents has necessitated the need for the rapid detection and identification of 4 common biological agents: *Bacillus anthracis*, *Yersinia pestis*, *Variola* and *Francisella tularensis*. Here we describe the utility of an integrated PCR-Microarray platform, the VereID™ Biosystem with the VereThreat™ Chip, which can enhance biological agent surveillance efforts. The platform is designed around a miniaturized silicon chip that integrates a PCR reactor with a customizable low-density DNA Microarray. As a result, the VereID™ Biosystem platform is able to perform multiplexing PCR and detection of the PCR products on a low-density DNA microarray, on the same chip. The probes of the microarray were designed to identify any or all 4 biological agents in one single test.

Materials and Methods


Positive controls preparation. The target genes of the biological agents are: a) *rpoB*, *pag* and *capB* genes for *Bacillus anthracis*; b) *pla* and *caf1* genes for *Yersinia pestis*; c) HA gene for *Variola virus* and; d) *tul4* and *fopA* genes for *Francisella tularensis*. The target region of the genes (~200bps) were synthesized, cloned into plasmids and used as surrogates for the analytical specificity and sensitivity tests.

PCR. For each biological agent, a PCR reaction of 25µl was setup with 10² copies of its corresponding genes and pipetted into the VereThreat™ Chip. The chip is placed into the VereID™ Biosystem for PCR amplification.

Chip Hybridization. Hybridization buffer was added to the amplicons and hybridization occurred for 30min at 62°C. After hybridization, chips were washed in SSC buffers and dried before imaging.

Microarray Imaging and Analysis. Hybridized samples were scanned using the Optical Reader. Signals were considered to be positive when its fluorescence intensity was more than three times that of the background.

VereThreat™ Chip Probe Layout



Capture probes for:

- *B. anthracis* ● *Y. pestis* ● *F. tularensis* ● *Variola virus*
- Orientation Probes ● Positive control ● Negative Control
- Hybridization control ● Binding Probes ● Probe location left empty

VereID™ VereThreat™ Chip operational procedure





FIG. 1 Typical microarray images demonstrating correct detection of (A) *B. anthracis*, (B) *Y. pestis*, (C) *F. tularensis* and (D) *Variola* virus. (E) Negative control without DNA.

Result and conclusion

Detection of biological agents. Analysis of the fluorescent images of the microarray showed that the VereThreat™ Chip was able to identify and discriminate the 4 targeted biological agents.

No cross-reactivity of the hybridization probes between the targets was observed.

With this Lab-on-Chip approach, results can be obtained rapidly, within 2 hours and is an invaluable tool for the rapid identification of *Bacillus anthracis*, *Yersinia pestis*, *Variola* and *Francisella tularensis*.

References

- Baldi, P., and G.W Hatfield. 2002. DNA microarrays and gene expression. From experiments to data analysis and modeling. University Press, Cambridge
- Keiko Tomioka, Michael Peredelchuk, Xiangyang Zhu, Roberto Arena, Dmitri Volokhov, Angamuthu Selvapandiyar, Katie Stabler, Jenny Mellquist-Riemenschneider, Vladimir Chizhikov, Gerardo Kaplan, Hira Nakhasi and Robert Duncan . 2005. A Multiplex Polymerase Chain Reaction Microarray Assay to Detect Bioterror Pathogens in Blood. J. of Molecular Diagnostics, Vol 7, No. 4: 486-94