

Foodborne Pathogen Detection and Identification Using VereFoodborne™ Chip Technology

Mitsuharu Sato, Dawn Liu and Rosemary Tan

Veredus Laboratories Pte Ltd,

83 Science Park Drive #04-02, The Curie, Singapore Science Park 1, Singapore 118258


Introduction

Foodborne diseases are responsible for a wide range of illnesses and are a growing public health problem all over the world. Current methodology for foodborne identification requires much time and effort and this slows down the response of health authorities.

Here, we describe the technology of an integrated PCR-microarray platform, the VereFoodborne™ Chip. This platform is based on a miniaturized silicon chip that integrates a PCR reactor and DNA microarray. Hybridization of the multiplex-PCR products by DNA microarray enables multiple, accurate and highly sensitive detection of foodborne pathogens in a single test.

#	pathogen	target gene		
1	Vibrio cholerae	Tu		
2	Vibrio parahaemolyticus	Tu		
3	Staphylococcus aureus	Tu		
4	Listeria monocytogenes	Tu		
5	Bacillus cereus	Tu		
6	Clostridium perfringens	Tu		
7	Campylobacter jejuni/coli/lari	Tu		
8	Escherichia coli	mutS	Stx 1A and 2A	STEC
9	Shigella flexneri/sonnei/dysenteriae	mutS+sp		S. dysenteriae
10	Salmonella	mutS	ttr, invA	
11	Cronobacter sakazakii	sp		
12	Norovirus	pol-capsid		
	PCR control	Luciferase		

Table 3 VereFoodborne™ target pathogens (ver. 3.0).
Tu: elongation factor Tu; mutS: mismatch repair protein sp: specific region found in our sequence analysis.



	Result	Sample information (showed after the test)
Test A	Positive for Vibrio parahaemolyticus, Campylobacter and Salmonella.	Vibrio parahaemolyticus (50 CFU eq.), Campylobacter coli (50 CFU eq.), Salmonella spp. (5 CFU eq.).
Test B	Positive for Vibrio parahaemolyticus, Campylobacter and Salmonella.	

Figure 5 VereFoodborne™ chip result 2 (mixed sample). Mixed genome sample was given from National Environment Agency, Singapore (NEA). The sample was tested on VereFoodborne version 1.1 and 2.0.



Figure 1 VereFoodborne™ assay process. 5 chips are processed in one TCS. Assay time (steps 1 to 6): 2 hrs




Figure 4 VereFoodborne™ chip result 1 (single target). 300-1,000 copies of genomic DNA or 100 copies of plasmid (stx1A) was tested as a template. Target signals are marked by color circles.




Figure 6 Experimental scheme of Salmonella spiking test. Salmonella was inoculated onto separated tomato and lettuce samples, then enriched with universal pre-enrichment broth as suggested in Bacteriological Analytical Manual (BAM) procedure. Salmonella detection was performed by BAM and VereFoodborne chip.

Test #1

Method	VereFoodborne	FDA BAM
Dosing bacteria	S. Enteritidis ATCC 13076 (Salmonella Group D)	
Dosing level	8.7 cfu/g	
Salmonella	mutS gene	34.8% (8/23*)
	invA gene	0% (0/23*)
	ttr gene	43.5% (10/23*)
Results	43.5% (10/23)	16.7% (4/24)
Detection time (hours)	2 - 3	72

The following bacteria were not spiked but were simultaneously detected in the VereFoodborne™ Chip

Bacteria	VereFoodborne	FDA BAM
E. coli	95.5% (21/23*)	Not tested
L. monocytogenes	47.8% (11/23*)	Not tested
B. cereus	4.3% (1/23*)	Not tested
C. perfringens	0% (0/23*)	Not tested

Test #2

Method	VereFoodborne	FDA BAM
Dosing bacteria	S. Typhimurium ATCC 14028 (Salmonella Group B)	
Dosing level	0.4 cfu/g	
Salmonella	mutS gene	63.6% (14/22)
	invA gene	86.4% (19/22)
	ttr gene	95.5% (21/22)
Results	95.5% (21/22)	100% (22/22)
Detection time (hours)	2 - 3	72

The following bacteria were not spiked but were simultaneously detected in the VereFoodborne™ Chip

Bacteria	VereFoodborne	FDA BAM
E. coli	81.8% (18/22)	Not tested
L. monocytogenes	13.6% (3/22)	Not tested
B. cereus	9.1% (2/22)	Not tested
C. perfringens	18.2% (4/22)	Not tested

Test #3

Method	VereFoodborne	FDA BAM
Dosing bacteria	S. Typhimurium (Salmonella Group B)	
Dosing level	0.4 cfu/g	
Salmonella	mutS gene	95.5% (21/22)
	invA gene	100% (22/22)
	ttr gene	100% (22/22)
Results	100% (22/22)	100% (22/22)
Detection time (hours)	2 - 3	72

The following bacteria were not spiked but were simultaneously detected in the VereFoodborne™ Chip

Bacteria	VereFoodborne	FDA BAM
E. coli	40.9% (9/22)	Not tested
L. monocytogenes	4.5% (1/22)	Not tested
B. cereus	9.1% (2/22)	Not tested
C. perfringens	4.5% (1/22)	Not tested

Figure 7 Salmonella spiking test results.

Test #1, #2 and #3 were conducted by National Environment Agency of Singapore (NEA) and Agri-Food and Veterinary Agency of Singapore (AVA), respectively.

Conclusions

In the Salmonella spiking test, the VereFoodborne chip showed similar sensitivity with the Bacteriological Analytical Manual (BAM) method in lettuce, and higher sensitivity in tomato samples. Highlighted differences between the two methods is *sample-to-answer* time and multiple target detection from one test. The assay time of VereFoodborne chip is 2-3 hours after enrichment, while that of BAM is 72 hours as well as its ability to detect multiple targets. This technology will significantly reduce the reaction time needed to address a potentially deadly threat.

Acknowledgements:

We would like to thank the Singapore National Environment Agency and the Singapore Agri-Food and Veterinary Authority for participating in this study.

VereFoodborne™ Panel 3.0

1	2	3	4	5	6	7	8	9	10	11
AT683	BA_BC1s	SHK_2td	ttr_2s	G1_157EC3h	BA_SA2s	G15TX_3s	BA_SA1s_S1D2	G2_VP2s	empty	AT809
AT730	G1_SE2s	LUC_2s	E_C1s	AT809	NLVG1_3s	STK2A_3	SAKA_4s	SHX_4sd	G15TX_2s	empty
empty	STK2A_1	BA_CP1s	G1_SHK1s	BA_BC3s	AT683	BA_LM4s	G1_EC5s	invA_4s	empty	empty
NLVG1_1s	LUC1_1s	NLVG1_2s	G2_VP1s_plus	SHX_3s	NLVG2_2s_pH	empty	empty	LUC3_3s	AT809	empty
AT809	NLVG2_1s	G1_1257EC2s	invA_2s	G1_EC4s	LUC1_2s	SHX_4s	LUC_3s	ttr_4s	AT776	empty
AT683	STK2A_2	G2_VC1s	BA_LM2s	SAKA_2s	BA_CP2s	G1_SE3s	E_C12s	G1_157SHX_3s	empty	AT809

12	13	14	15	16	17	18	19	20	21
empty	BA_BC1s	SHK_2td	ttr_2s	G1_157EC3h	BA_SA2s	G15TX_3s	BA_SA1s_S1D2	G2_VP2s	AT683
AT730	G1_SE2s	LUC_2s	E_C1s	AT809	NLVG1_3s	STK2A_3	SAKA_4s	SHX_4sd	G15TX_2s
empty	STK2A_1	BA_CP1s	G1_SHK1s	BA_BC3s	AT683	BA_LM4s	G1_EC5s	invA_4s	empty
NLVG1_1s	LUC1_1s	NLVG1_2s	G2_VP1s_plus	SHX_3s	NLVG2_2s_pH	empty	empty	LUC3_3s	AT809
AT809	NLVG2_1s	G1_1257EC2s	invA_2s	G1_EC4s	LUC1_2s	SHX_4s	LUC_3s	ttr_4s	AT776
empty	STK2A_2	G2_VC1s	BA_LM2s	SAKA_2s	BA_CP2s	G1_SE3s	E_C12s	G1_157SHX_3s	AT683

Figure 2 VereFoodborne™ probe layout (ver. 3.0).

Layout is in 6 x 21 format, which consists of:
84 capture probes (42 x 2 replicates)
10 hybridization control probes
10 PCR control probes
8 hybridization negative control probes