

## Protective Face Shield Visor – Reference 2002



Our Ref 2002 is a lightweight, clear face shield for personal use, that will create a barrier to help reduce droplets and sneezes.

### USP

- The Protective Face Shield Visor – Reference 2002 is lightweight in use and easy to wear
- The Face Shield can be cleaned with warm water and any kind of cleaning detergent

T: (0191) 269 7810 | F: (0191) 267 1028

E: [sales@ukindustrialtapes.co.uk](mailto:sales@ukindustrialtapes.co.uk) | [sales@postapes.co.uk](mailto:sales@postapes.co.uk) | [sales@stix2.co.uk](mailto:sales@stix2.co.uk)

Brumwell House, Westway Industrial Estate, Throckley, Newcastle Upon Tyne, NE15 9EW

[www.ukindustrialtapes.co.uk](http://www.ukindustrialtapes.co.uk) | [www.postapes.co.uk](http://www.postapes.co.uk) | [www.stix2.co.uk](http://www.stix2.co.uk)



It is the buyer's responsibility to consider & assess the fitness for purpose of the goods, having regard to use to which the buyer intends to put such goods. The price is quoted for the quantity quoted. All prices quoted are for 30 days only. Standard credit terms are 30 days from the end of month unless agreed in writing otherwise. All orders are subject to this Company's terms & conditions of trading which can be seen on our website [ukindustrialtapes.co.uk](http://ukindustrialtapes.co.uk), [postapes.co.uk](http://postapes.co.uk) and [stix2.co.uk](http://stix2.co.uk)



POS Tapes & Stix2 are trading divisions to UK Industrial Tapes Ltd

### Application

- Personal use
- Help support essential business and NHS staff

### Product Description

- A lightweight, clear face shield for personal use, that will create a barrier to help reduce droplets and sneezes. The parts used are removable, cleanable and replaceable.

### Features

- Very lightweight material
- Practical minimal design
- Reduction of droplets from coughs and sneezes

### Application

- Personal use
- Help support essential business staff

### Advantages

- Extremely lightweight material
- Easily cleaned
- Reusable
- Long service life use

### Storage Conditions

- When stored in stacks, protect against heat, moisture and direct sunlight
- Do not expose the Face Shield with the protective film to the outside weather

### Use & Handling

- Instructions on how to prepare and use the Face Shield have been provided
- Place your hand on the inside when putting on the Face Shield to avoid getting scratched by the edges

### Cautions

- Do not apply much pressure when bending the Face Shield to avoid breaking it

### Packaging

They are packed separately in packs of 10.  
10 straps and 10 visors in a sealed bag  
Master boxes of 250 visors complete with unbranded instructions.

### Product Material

- Protective Face Shield Visor – Reference 2002 (456.1 x 240)
- 275 clear Micron APET

\*Disclaimer: Whilst our Protective Face Shield Visor - Reference 2002 has many applications it is not designed to be used with Covid -19 positive patients. However, we are currently working on a successor to our 2002 which is currently being tested in the clinical setting and CE approval has already been applied for.

T: (0191) 269 7810 | F: (0191) 267 1028

E: [sales@ukindustrialtapes.co.uk](mailto:sales@ukindustrialtapes.co.uk) | [sales@postapes.co.uk](mailto:sales@postapes.co.uk) | [sales@stix2.co.uk](mailto:sales@stix2.co.uk)

Brumwell House, Westway Industrial Estate, Throckley, Newcastle Upon Tyne, NE15 9EW

[www.ukindustrialtapes.co.uk](http://www.ukindustrialtapes.co.uk) | [www.postapes.co.uk](http://www.postapes.co.uk) | [www.stix2.co.uk](http://www.stix2.co.uk)

