

# Fevin Scan Skin Temperature Measurement System





MIRICLE FevIR Scan 110K Thermal



Multi Target Auto Tracking System



High Accuracy Infrared Thermal Calibration Source



### Fev/R<sup>™</sup> Scan Skin Temperature Measurement System

Ideal for Mass Screening. In today's world of fast transit global travel, outbreaks of infectious disease will cross national and international boundaries in a matter of hours. The SARS, Bird 'Flu and Swine 'Flu,  $H_1N_1$  and now Ebola virus, have affected tens of thousands of people and businesses with massive disruption to all people in all walks of life.

Miniature, discreet and ultra portable MIRICLE FevIR Scan Skin Temperature Monitoring system provides peace of mind with high speed sensitive temperature measurement for high volume pedestrian transit areas.

MIRICLE FevIR Scan thermal imaging and temperature measurement system from Thermoteknix is easily installed, operated and re-located with the very minimum of setup time and operator training.

On screen visible and audible alerts provide instant notification of high temperature alarms. On-screen alerts identify one or more individuals in a large crowd situation, aiding rapid dectection and maintaining passenger flow.

The Thermoteknix MIRICLE 110K radiometric infrared thermal imaging camera can detect temperature differences as small as 0.05°C and with its high resolution 110,000 pixel detector and lens options, covers a wide or narrow field of view as required for transit monitoring.

Password protected configuration to prevent unauthorised changes.

Operating at 50Hz (9Hz licence free models available) MIRICLE 110K ensures that a high rate of human transit is monitored effectively without slowing queues in busy areas.

## MIRICLE FevIR Scan Skin Temperature Measurement System: Accurate, Affordable - Quick and Easy to Install - Peace of Mind


Caution: Human Skin temperature varies with environment and physiological factors and does not always accurately reflect internal core temperature. Raised Skin temperatures can occur in the absence of fever and vice versa.

#### **European Headquarters**

Teknix House, 2 Pembroke Avenue Waterbeach, Cambridge, CB25 9QR UK Tel: +44 (0)1223 204000 Fax: +44 (0)1223 204010 sales@thermoteknix.com www.thermoteknix.com



# Fevir Scan Skin Temperature Measurement System













## MIRICLE FevIR Scan Skin Temperature Measurement System Specification

MODEL	110KR
Detector material	Amorphous Silicon (A Si)
Array size	384x288
TEC-Less operation	YES
Pixel count	110,592
Spectral response	8µm-12µm
Pixel Pitch	25µm
Response time	<10 mS
Sensitivity (Excluding optics)	0.05°C
Frame rate	50/60 Hz (9 Hz option available)
Video Output	USB 2
Power consumption	2.5 Watts
Time to image	≤2 Seconds
Std. ambint Operating temperature range	-20 to + 50°C
Remote command interface	RS232
PC Interface	USB 2 - 50 Hz

CONFIGURATION	MIRICLE 110KR
Weight (Excluding optics)	166g
Length (mm)	57.6
Width (mm)	45
Height (mm)	52.5

FevIR Scan is available without export licence with a 9Hz detector update rate. 50/60Hz update rate cameras may be subject to export licence for some countries

Installation Options:

- 1. Tripod Pack
- 2. Frame Pack





Caution: Human Skin temperature varies with environment and physiological factors and does not always accurately reflect internal core temperature. Raised Skin temperatures can occur in the absence of fever and vice versa. Specification Subject to change without notice.

**European Headquarters**, Teknix House, 2 Pembroke Avenue Waterbeach, Cambridge, CB25 9QR UK Tel: +44 (0)1223 204000 Fax: +44 (0)1223 204010 sales@thermoteknix.com www.thermoteknix.com