PowerLaser PRO[™]

Probably the best handheld laser in the world



Accelerated Natural Healing[™]

Biostimulating laser therapy is a unique and unmatched modality working on the body's own terms with no known side effects.

Laser therapy stimulates the body's own natural repair system and accelerates all phases of the natural healing process by:

- · Enhancing the entire immune system
- · Stimulating cell proliferation
- · Increasing blood circulation
- · Increasing formation of collagen tissue
- Relieving pain
- · Increasing the production of ATP



With acute injuries, like sprains and other sports injuries, healing time is reduced dramatically. Full recovery is reached within a fraction of the time it would normally take, left to natural healing. Healing time is typically reduced "from weeks to days".


Applications of laser therapy

Physiotherapy Chiropractic Dentistry Veterinary General practice Wound care Dermatology Beauty treatment Acupuncture


HyperPulsed[™]

The PowerLaser PRO Series comes with a brand new, proprietary technology that we call HyperPulsed[™].

The energy is delivered in high peak pulses of a staggering 1000 mW per diode, pulsed with a duty cycle of 50%, so the resulting average output power is 500 mW per diode.



This will give a much deeper penetration into the tissue, compared to 500 mW continuous wave, while maintaining an average output power that prevents excessive heating of the tissue.



Despite the deep penetrating peak power output of 1000 mW, the average output of 500 mW keeps the laser in Class 3B.

It is a Class 3B laser with the punch of a Class 4 laser, lightweight, portable and totally independent.

Indications

PowerLaser PRO[™] is CE certified for: Arthritis Inflammation Epicondylitis Myofascial pain Wound healing

PowerLaser PRO[™] is FDA cleared for pain relief, carpal tunnel syndrome and wound healing.

The PowerLaser PRO[™] advantages

Touch control

PowerLaser PRO[™] is controlled by your fingertips - just slide along the touch panel for program select and start / stop of the laser.

9 different protocols

PowerLaser PRO[™] has 3 different output power settings and 3 different exposure time settings, giving a total of 9 different protocols. The combinations are carefully selected according to the latest scientific research, for precise adjustment to each individual therapy situation.

Wireless

PowerLaser PRO[™] is battery operated and independent of power supply or wall outlet.

Heavy duty design

Sturdy and durable, all metal Scandinavian design, withstands wear and tear of daily use.

Vibration

A discrete vibration every 2 sec. when the laser is active provides tactile feedback to practitioner and patient. You can **feel** that the laser is ON.

High quality sound

Each program choice has its own, individual tone. The audible signal is round and pleasant. You can **hear** what your laser is doing.

The perfect laser therapy hand tool

Professional practitioners worldwide praise two distinct features of the PowerLaser PRO[™]:

- 1. the portability
- 2. the ease of use

"I carry the PowerLaser in my pocket all the time at the clinic. I use it on nearly all my patients. It is adjusted and ready to go within seconds. No dialing through endless menus. Only 3 power levels, from the superficial to the very profound. It can't be any easier."

Approvals

The PowerLaser PRO[™] Series is produced under the highest industrial standards in a ISO9001 certified environment.

The production conforms to following standards and specifications:

Annex V, section 3.2 - Production Quality Assurance. EN 60601-2-22:1996.93/42/EEC concerning medical devices. Class 2a medical devices. Class 3B laser products. The PowerLaser PRO is available in 2 models: PLp1500: 1500mW HyperPulsed to 3000mW peak PLp500: 500mW HyperPulsed to 1000mW peak

Ask for pricing and detailed information:

Quantum-Healing-Lasers.com 1215 S Kihei Rd Ste O #810 Kihei, HI 96753 Serving all of North America

- T: 888-824-7558
- E: contactus@coldlasersupplies.com
- W: Quantum-Healing-Lasers.com