



Cavendish House Dental Practice | Case Study



Cavendish Dental Practice were looking to implement the use of a new patient camera and DVD player into one of their surgery rooms. However, there was no obvious space in which to store the new system. Restrictions prevented the monitor from being mounted off the wall or floor. Additionally to this, the new system would be required to adjust to a comfortable viewing position to suit both the patient and the dentist.

Internal cable management was also key in order to reduce dust collection points. After careful assessment of the dental theatre, it was decided that a height and reach adjustable ceiling mounted solution would provide all the movement required. This space saving ceiling mount solution also has the added benefit of zero floor footprint.

After defining key objectives within the Cavendish Dental Practice it was decided that the T2 Elite Ceiling Mount would exceed all criteria's. The EMT2EQ-C812- QG provides a massive 1447mm of reach and 533mm of gas assisted height adjustment. All cables are routed internally through each arm. This solution allows both dentist and patient to comfortably view the monitor and adjust between positions with ease. ErgoMounts completed the installation on time and on budget.

Part List:

EMT2EQ-C8X12-QG



T2 Elite Ceiling Mount 12″ Drop Light Grey

Additional Services:



T +44 (0)1252 333326 **F** +44 (0)1252 330006 **E** sales@ergomounts.co.uk

see it your way...

www.ergomounts.co.uk