

### STERILE LINERS

- Easy to install
- Complete coverage of tray top and sides:
  - liners attach at tray bottom
- Liners clearly identify **Clean** vs **Biohazardous** 360° visibility™
- Compatible with other transport trays, bins, and containers

*SafeCAP saves us the cost of torn liners, positive ATP tests and additional cleaning.*

GI Lab Manager  
300 bed hospital  
Southern California

*SafeCAP is easier to install and provides better patient and staff protection.*

GI Director  
Midwest Healthcare System



 **CLEAN**



 **BIOHAZARDOUS**

### TRANSPORT TRAY AND LIDS

- Reusable, washable at high temperatures
- Elegant design: no difficult to clean areas
- Ergonomic shelf for easy handling
- Safe transport
- Solid, secure lid allows for stacking of trays
- Compatible with other transport carts
- Meets SGNA®, AORN® and AAMI® Guidelines
- Meets OSHA® requirements

Easy to clean!



# SafeCAP®

by ClinicalChoice

Meets SGNA®, AORN®, AAMI®, and OSHA® Guidelines

**SafeCAP®** is a user-friendly endoscope transport and short-term storage system. Easy to install Clear, green **Clean** and red **Biohazardous** liners are compatible with other transport trays, containers and bins. Easy to clean and disinfect trays and lids are compatible with most carts and cart washers. **SafeCAP®** makes it simple to practice better infection control and experience immediate cost savings.



## SafeCAP® Transport Cart

- Easy access
- Locking wheels
- Clear visibility
- Secure tray transport

Features	✓ SafeCAP®	✗ Competitive System	✗ Bag System	✗ Disposable Tray System
😊 <b>CLEAN</b>   ☠️ <b>BIOHAZARDOUS</b> 360° visibility™	✓	✗	✗	✗
Cart washer compatible trays and lids*	✓	✗	✗	✗
Solid lid design allows for safe stacking of trays	✓	✗	✗	✗
Liners cover trays sides and attach at the bottom	✓	✗	✗	✗
Contamination kept in/out on top and sides	✓	✗	✗	✗
Easily read 😊 <b>CLEAN</b> and ☠️ <b>BIOHAZARDOUS</b> liners on carts	✓	✗	✗	✗
Easy to clean and disinfect tray and lid	✓	✗	✗	✗
Tray supports weight of endoscope	✓	✓	✗	✗
Individually packaged liners	✓	✗	✗	✗
Eliminates areas for bacteria to hide	✓	✗	✗	✗

ORDERING INFORMATION		
Part Number	Product Description	Units Per Box
SC-1478	Sterile Liner Set (Clear, 😊 <b>CLEAN</b> and ☠️ <b>BIOHAZARDOUS</b> liners)	100
SC-CLR-100	Clear Sterile Liners - 100 each individually packaged	100
SC-GRN-100	😊 <b>CLEAN</b> Sterile Liners - 100 each individually packaged	100
SC-RED-100	☠️ <b>BIOHAZARDOUS</b> Sterile Liners - 100 each individually packaged	100
SC-152L	Clear Reusable Lid**	1
SC-249G	Gray High Temperature Tray*	1
SC-249L	Gray High Temperature Lid*	1
SC-CART	Seven tray cart**	1

\*Maximum temp. 230°F/110°C

\*\*Maximum temp. 180°F/82°C

### SafeCAP® Endoscope Transport and Short-Term Storage System | © 2019 and beyond

The following patents are associated with SafeCAP® products: U.S. Pat. No. D818,841 | U.S. Pat. No. D819,409

SGNA®, AORN®, AAMI® and OSHA® are registered trademarks of Society of Gastroenterology Nurses Association and Associates, Inc., Association of Operating Room Nurses, Inc., Association for Advancement of Medical Instrumentation and U.S. Department of Labor, Occupational Safety and Health Administration.