

LUMIPULSE® G1200



SIMPLIFY YOUR WORK USING A FULLY AUTOMATED IMMUNOASSAY SOLUTION



**A robust market proven CLEIA analyzer
utilising mono-test cartridges**

LUMIPULSE is a registered trademark of Fujirebio Inc.



LUMIPULSE® G1200

FULLY AUTOMATED CHEMILUMINESCENT ENZYME IMMUNOASSAY SYSTEM

01 CONSTANT THROUGHPUT

- 120 tests/hour independent of assay combinations or assay format order
- Each 30 seconds a result
- Automated dilution



02 CARTRIDGE

- Interchangeable on all LUMIPULSE G series
- Cost-transparent
- No wastage from on-board expiry
- Simple operation



03 USER INTERFACE

- Simple menu navigation
- Limited hands-on time



PROVEN RELIABILITY


Meantime between failure
1125 days (figures April 2017)



04 NO SPECIMEN CARRY-OVER
 Instrument uses disposable sample tips equipped for level sensing



05 INCREASED WALK AWAY TIME
 Load up to 504 tests and 100 samples



06 STABLE, SIMPLE CALIBRATORS
 Most assays use 2-3 point calibrations and stable for 30 days



READY TO CONNECT TO LABORATORY AUTOMATION SYSTEMS

System architecture optimized for connection to automated sample tracking systems

LUMIPULSE TRACK CONNECTION



© Fujirebio Europe, May 2017, FRE-027, r1

*For designing your
personal solution,
please contact
your local
representative.*

Fujirebio Europe N.V.
Technologiepark 6 - 9052 Gent - Belgium
T +32 (0)9 329 13 29 - F +32 (0)9 329 17 75
www.fujirebio-europe.com

