

Price List

Valid as of 2015-09-01

The volatile exchange rates of the euro and the US dollar have made it impossible to have fixed prices in these currencies. We will, however, try to maintain fixed prices in Swedish Crowns (SEK), which actually have been fixed since 2012-04-01. You may order and pay in EUR or USD, but the price will be based on the exchange rate, which will be updated on a regular basis. As a service to our customers indicative prices in EUR and USD are shown in the two right-hand columns. Actual prices are based on the prices in SEK. This may lead to slight discrepancies between the price at the time of order and that at the time of order confirmation. We ask for your understanding in this matter.

Art. No.	Description	Notes	Package	Conc.	Size	Price SEK	Indicative price		
							EUR	USD	
							Exchange rates versus SEK		
							11,09	10,37	
Kits									
111-021	Kinase RR Kit		kit		1000 tests	2596	234	250	
144-041	ATP Kit SL		kit		1000 tests	3207	289	309	
155-050	Cellular ATP Kit HTS		kit		1200 tests	3968	358	383	
188-441	Cell Viability Kit SL		kit		1600 tests	6943	626	670	
266-111	Intracellular ATP Kit HS		kit		375 tests	4125	372	398	
266-112	Microbial ATP Kit HS		kit		375 tests	4218	380	407	
266-221	ATP Hygiene Kit HS (dropper bottles)	a, b	kit		100 tests	2039	184	197	
266-311	ATP Biomass Kit HS		kit		375 tests	3365	303	324	
484-001	Luciferase Assay Kit		kit		800 tests	3662	330	353	
484-102	Luciferase Assay Kit HTS		kit		4800 tests	18364	1656	1771	
484-103	Luciferase Assay Kit HTS	c	kit		4600 tests	18364	1656	1771	
PHL-192	Portable Hygiene Laboratory	d	kit		192 tests	3699	334	357	
Reagents									
11-501	ATP Reagent SL		vial		250 tests	612	55	59	
11-501-M	ATP Reagent SL (excl. Mg2+)		vial		250 tests	612	55	59	
12-101	ATP Reagent HS		vial		100 tests	343	31	33	
12-071	ATP Reagent HS		vial		17 tests	260	23	25	
13-101	ATP Reagent SS	a, e	vial		100 tests	1835	165	177	
51-101	ATP Substrate		vial		100 tests	232	21	22	
51-201	ATP Substrate		vial		200 tests	426	38	41	
61-101	Luciferin Substrate		vial		100 tests	232	21	22	
61-201	Luciferin Substrate		vial		200 tests	426	38	41	
71-101	ATP Eliminating Reagent		vial		500 tests	426	38	41	
BT11-100	D-Luciferin		vial		100 mg	458	41	44	
BT11-500	D-Luciferin		vial		500 mg	1465	132	141	
BT11-1000	D-Luciferin		vial		1000 mg	2345	211	226	
BT12-100	L-Luciferin		vial		100 mg	1608	145	155	
BT11-100 K	D-Luciferin K salt		vial		100 mg	458	41	44	
BT11-1000 K	D-Luciferin K salt		vial		1000 mg	2345	211	226	
BT11-100 Na	D-Luciferin Na salt		vial		100 mg	458	41	44	
BT11-1000 Na	D-Luciferin Na salt		vial		1000 mg	2345	211	226	
LA-10	Firefly Luciferase (recomb., enhanced activity)		vial		10 mg	7323	660	706	
For other luciferase products see under Kikkoman									
Standards									
47-051	ATP Standard		vial	100 nmol/L	5 mL	185	17	18	
46-051	ATP Standard		vial	1 µmol/L	5 mL	185	17	18	
45-051	ATP Standard		vial	10 µmol/L	5 mL	185	17	18	
48-072	ATP Standard (dropper bottle)		bottle	10 nmol/L	6.7 mL	278	25	27	
Buffers, Diluents, etc									
21-101	Tris-EDTA Buffer		vial	0.1 mol/L	10 ml	83	7	8	
21-103	Tris-EDTA Buffer		vial	0.1 mol/L	100 ml	417	38	40	
21-500	Tris-EDTA Buffer		bottle	0.1 mol/L	500 ml	1353	122	130	
22-101	Diluent B		vial		10 ml	185	17	18	
22-072	Diluent B (dropper bottle)		bottle		6.7 ml	158	14	15	
23-101	Diluent C (SL)		vial		10 ml	185	17	18	
24-101	Diluent D		vial		10 ml	185	17	18	
27-101	ATP Free H2O		vial		10 ml	176	16	17	
40-500	ATP Diluent		vial		10 ml	222	20	21	
Extractants									
31-072	Extractant B/S (dropper bottle)		vial		6.7 ml	213	19	21	
31-101	Extractant B/S		vial		10 ml	315	28	30	
31-103	Extractant B/S		bottle		100 ml	1613	145	156	
31-500	Extractant B/S		bottle		500 ml	4867	439	469	
20-121	Lysing Diluent		vial		12 ml	213	19	21	
28-101	Cell Lysing Reagent		vial		10 ml	473	43	46	
Disposables									
681000	Cuvettes with 'stirring wings'		bag		1000 pcs	1261	114	122	

Kikkoman (for European customers only)**Reagents**

61303	LuciPac W	f	kit	100 tests			
60331	LuciPac Pen	f	kit	100 tests			
60336	LuciPac Pen- Agua	f	kit	100 tests			
60361	LuciPac A3 Surface	f	kit	100 tests			
60365	LuciPac A3 Water	f	kit	100 tests			
	CheckLite 250 Plus		kit	250 tests	3454	311	333
	CheckLite HS Set		kit	100 tests	3454	311	333
	CheckLite CT150		kit	150 tests	4001	361	386
LA-11	Firefly Luciferase (recomb., thermostable)		vial	10 mg	3232	291	312
LA-12	Firefly Luciferase (recomb., thermostable)		vial	100 mg	13150	1186	1268

Systems

PD-20	Lumitester™ handluminometer PD-20	f					
PD-30	Lumitester™ handluminometer PD-30	f					
61234	Lumitester Smart	f					
C-110	Lumitester™ portable luminometer C-110 inc	f					
	1. LUMITESTER C-110						
	2. AC adapter						
	3. Carrying case (Cloth type)						
	4. Cable for PC						
	5. English version of 'Job Manager'						
DPU-H245	Printer set for LUMITESTER C-110	f					
	Includes following items						
	1. Printer for LUMITESTER C-110						
	2. AC adaptor						

Prices do not include VAT, shipping, customs and other duties and fees.

In almost all kits and reagents, the number of tests per kit is calculated based on using all the reagent (no dead volumes, etc.) and microplate use. The use of cuvettes yields 4 or 5 times less tests per kit depending on the product. This also applies to standard and blanks. Also see Notes below.

Notes:

- a For this product, just an application for cuvettes is available and the number of tests/kit is calculated based on pipetting accessing all reagents. Standard and blanks are regarded as tests.
- b Does not include swabs nor cuvettes
- c For dispensing luminometers.
- d Contains 266-113 ATP Hygiene Reagents (similar to 266-221 ATP Hygiene Kit HS) but also swabs, small cuvettes (microtubes), tips and even a pipette. The only item not included is a luminometer. The swabs are flat making swabbing much easier and more reliable. Intended for testing of medical devices, food contact surfaces and general cleaning control.
- e The most sensitive ATP Reagent on the market.
- f Subject to fluctuations of exchange rate JPY/SEK

BioThema AB
 Handens stationsväg 17
 SE-136 40 Handen, Sweden
 tel + 46 (0)8 7770020
www.biothema.com